

Target™ Detachable Coils

Target 360 Detachable Coils

Conformability from frame to finish


Open loop configuration designed to conform to different aneurysm shapes, fill uniformly, and minimize compartmentalization.



Target 3D Detachable Coils

Consistent framing

Alternating large and small loops are designed to create a stable frame.



Target XL Detachable Coils

Larger and longer

Target XL Coils provide two times more fill and 40% more width than Target 10-size coils. All deliverable through Excelsior SL-10 and XT-17 Microcatheters.



Target Nano Detachable Coils

Beyond soft

2x softer than Target Ultra Coils with advanced microcatheter stability. Designed to fill small spaces and conform to the final spaces within aneurysms.



InZone Detachment System

Smarter detachment

Electrolytic detachment affords a small and flexible coil junction, enhancing coil delivery.



Target Detachable Coil

See package insert for complete indications, contraindications, warnings and instructions for use.

Intended use / indications for use

Target Detachable Coils are intended to endovascularly obstruct or occlude blood flow in vascular abnormalities of the neurovascular and peripheral vessels.

Target Detachable Coils are indicated for endovascular embolization of:

- Intracranial aneurysms
- Other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae
- Arterial and venous embolizations in the peripheral vasculature

InZone Detachment System

See package insert for complete indications, contraindications, warnings and instructions for use.

Indications for use

The InZone Detachment System is intended for use with all versions of Stryker Neurovascular detachable coils in the embolization of intracranial aneurysms and other vascular malformations of the neuro and peripheral vasculature.

This document is intended solely for the use of healthcare professionals.

A physician must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that physicians be trained in the use of any particular product before using it in a procedure. The information presented is intended to demonstrate the breadth of Stryker product offerings. A physician must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

stryker

Target™ Detachable Coils



Smooth and stable

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: InZone, Stryker, Target, Target Nano, Target XL. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.


Copyright © 2022 Stryker
AF001968 v3.0

NZ **New Zealand
Sponsor Address**
Stryker New Zealand Ltd
PO Box 17136
Greenlane 1546
Auckland
New Zealand

AUS **Australian
Sponsor Address**
Stryker Australia Pty Ltd
8 Herbert Street
St Leonards, NSW 2065
Australia


Stryker Neurovascular
47900 Bayside Parkway
Fremont, CA 94538
stryker.com
Date of Release: JAN/2022
EX_EN_IL

Selection guide



360 shape

Target XL 360 Standard Coils			Target 360 Standard Coils			Target 360 Soft Coils			Target 360 Nano Coils		
Secondary diameter (mm)	Length in introducer (cm)	Product number	Secondary diameter (mm)	Length in introducer (cm)	Product number	Secondary diameter (mm)	Length in introducer (cm)	Product number	Secondary diameter (mm)	Length in introducer (cm)	Product number
6	20	M0036106200	3	6	M0035463060	2	4	M0035472040	1	2	M0035421020
7	20	M0036107200		8	M0035463080	3	6	M0035473060	3	3	M0035421030
8	30	M0036108300		8	M0035464080		8	M0035473080	2	2	M0035421520
9	30	M0036109300	10	M0035464100	10		M0035473100	1.5	3	M0035421530	
10	40	M0036101040	4	15	M0035464150	3.5	10	M0035473510	4	4	M0035421540
12	45	M0036101240		15	M0035465150	4	6	M0035474060	2	3	M0035442030
14	50	M0036101450		20	M0035465200		8	M0035474080	2.5	4	M0035442540
16	50	M0036101650	15	M0035466150	10		M0035474100	3	4	M0035443040	
18	50	M0036101850	6	20	M0035466200	15	M0035474150	6	6	M0035443060	
20	50	M0036102050		30	M0035466300	4.5	12	M0035474510	3.5	6	M0035443560
22	50	M0036102250		15	M0035467150	10	M0035475100	5	10	M0035475100	
24	50	M0036102450	20	M0035467200	15	M0035475150	2		3	M0035422030	
			30	M0035467300	20	M0035475200	4		4	M0035422040	
			8	20	M0035468200	10	M0035476100	6	6	M0035422060	
				30	M0035468300	15	M0035476150	2.5	4	M0035422540	
				30	M0035469200	20	M0035476200	4	4	M0035423040	
			9	30	M0035469300	30	M0035476300	6	6	M0035423060	
				10	M0035461030	15	M0035477150	8	8	M0035423080	
				11	M0035461130	20	M0035477200	10	10	M0035423100	
			12	30	M0035461230	7	20	M0035477300	3.5	8	M0035423580
				13	M0035461330	30	M0035478200	4	6	M0035424060	
				14	M0035461430	30	M0035478300	3	8	M0035424080	
			15	30	M0035461540	8	20	M0035478200	4	8	M0035424100
				20	M0035461540	30	M0035479200	30	10	M0035424150	
				30	M0035461540	30	M0035479300	9	15	M0035424150	
			16	50	M0036121450	10	30	M0035471030	4.5	10	M0035424510
				50	M0036121650	11	30	M0035471130	5	10	M0035425100
				50	M0036121850	12	30	M0035471230	5	15	M0035425150
			20	50	M0036122050	13	30	M0035471330			
				50	M0036122050	14	30	M0035471430			
				50	M0036122050	14	30	M0035471430			



3D shape

Target 3D Detachable Coils		
Secondary diameter (mm)	Length in introducer (cm)	Product number
3	6	M0036153060
4	8	M0036154080
5	10	M0036155100
6	15	M0036156150
7	15	M0036157150
8	25	M0036158250
10	30	M0036151030



Helical shape


Target XL Helical Coils			Target Helical Nano Coils			
Secondary diameter (mm)	Length in introducer (cm)	Product number	Secondary diameter (mm)	Length in introducer (cm)	Product number	
5	15†	M0036115150	1	2	M0035431010	
6	20†	M0036116200		2	M0035431020	
7	20†	M0036117200		3	M0035431030	
8	30†	M0036118300	1.5	1	M0035431510	
9	30†	M0036119300		2	M0035431520	
10	40†	M0036111040		3	M0035431530	
			4	M0035431540		
			2	2	M0035452020	
				3	M0035452030	
				4	M0035452040	
			2.5	6	M0035452060	
				3	M0035452530	
				4	M0035452540	
			3	6	M0035452560	
				4	M0035453040	
				6	M0035453060	
			3.5	8	M0035453080	
Target Helical Ultra Coils			Secondary diameter (mm)	Length in introducer (cm)	Product number	
			2	1	M0035432010	
				2	M0035432020	
				3	M0035432030	
			2.5	4	M0035432040	
				6	M0035432060	
				8	M0035432080	
			3	3	M0035432530	
				4	M0035432540	
				6	M0035432560	
			3.5	4	M0035433040	
				6	M0035433060	
				8	M0035433080	
			4	10	M0035433100	
				6	M0035433560	
				8	M0035434060	
			4.5	8	M0035434080	
				15	M0035434150	
			4.5	8	M0035434580	

Target Coil Accessories	
Description	Product number
InZone Detachment System	M00345100950

Target™ Detachable Coils

Smooth and stable

Consistently smooth deployment and enhanced microcatheter stability.



† Comparable to Soft Coils.
Not all products are available for sale in all countries.